
















NONABSORBABLE SURGICAL SUTURES

LINEX*	Monofilament Polyamide
Virtues	: Being monofilament, Linex is remarkably smooth and gives excellent knot security. Possesses good tensile strength. Easy removal with no tissue adherence.
Material	: Polyamide 6 suture (Nylon). Blue/Black coloured for easy identification in the operating field.
Tissue Reaction	: Linex Suture is inert like stainless steel. Maintains tensile strength indefinitely in tissues. Does not support infection.
Applications	: General closure, Skin and Plastic Surgery. Also used for Ophthalmic and Micro Surgery.
Choice	: Loop, Needled and pre-cut sterile sutures conforming to USP standards.
Caution	: Not recommended for attachment of artificial prostheses in Cardiovascular surgery.



NEEDLE LENGTH	NEEDLE DESCRIPTION	SUTURE LENGTH	SUTURE SIZE IN USP						
			10-0	6-0(0.7)	5-0(1)	4-0(1.5)	3-0(2)	2-0(3)	
 Double needle	3/8 Circle spatulated micro point	38 cm	FS 3719						
	3/8 Circle Reverse Cutting	70 cm		FS 3320					
	3/8 Circle Reverse Cutting	70 cm			FS 3317				
	3/8 Circle Cutting	35 cm			FS 3315				
		70 cm				FS 3318	FS 3321		
	3/8 Circle Reverse Cutting	100 cm					FS 3328	FS 3329	
		45 cm					FS 3328S	FS 3329S	

NEEDLE LENGTH	NEEDLE DESCRIPTION	SUTURE LENGTH	SUTURE SIZE IN USP		
			2-0(3)	1-0(3.5)	1(4)
	3/8 Circle Cutting	100 cm	FS 3236	FS 3239	
	1/2 Circle Reverse cutting	100 cm		FS 3337	
		45 cm		FS 3337S	
	1/2 Circle Round Bodied	150 cm (Loop)		FS 3340	
	1/2 Circle Reverse cutting	100 cm		FS 3339	
	1/2 Circle Round Bodied	100 cm		FS 3346	FS 3347

NEEDLE LENGTH	NEEDLE DESCRIPTION	SUTURE LENGTH	SUTURE SIZE IN USP		
			2-0(3)	1-0(3.5)	1(4)
 45 mm	3/8 Circle Reverse Cutting	45 cm	FS 3336S		
		100 cm	FS 3336	FS 3335	FS 3334
 45 mm Heavy	1/2 Circle Reverse Cutting	100 cm			FS 3338
 50 mm Heavy	1/2 Circle Round Bodied	150 cm (Loop)			FS 3348
 60 mm	Straight Cutting	70 cm	FS 3390	FS 3389	
 90 mm	3/8 Circle Cutting	100 cm			FS 3397

PRE-CUT SUTURES

SUTURE LENGTH	PACKAGING	SUTURE SIZE IN USP		
		2-0(3)	1-0(3.5)	1(4)
5 x 35 cm	Sterilised Monofilament in sterilised foil over wrap. 5 x 35 cm box of 12's pack	FS 703	FS 704	FS 705
2 x 76 cm	Sterilised Monofilament in sterilised foil over wrap. 2 x 76 cm box of 12's pack	FS 903	FS 904	FS 905