




















NONABSORBABLE SURGICAL SUTURES

SUTURA*	Black Braided Silk
Virtues	: Most widely used. Excellent handling characteristics. Good knot security.
Material	: Entwined thread from silk worm larva. Treated with silicone or wax to give strength and smooth surface.
Tissue Reaction	: Moderate
Applications	: General closure, G. I. Tract, Plastic, Skin, Ophthalmic, Cardiovascular Surgeries. Widely used as ligature.
Choice	: Needled and pre-cut sterile sutures conforming to USP standards.
Caution	: Avoid using in the presence of infection. Not used in biliary and urinary tract surgery.



NEEDLE LENGTH	NEEDLE DESCRIPTION	SUTURE LENGTH	SUTURE SIZE IN USP			
			5-0(1)	4-0(1.5)	3-0(2)	2-0(3)
 12 mm	3/8 Circle Reverse Cutting	76 cm	FS 5080			
 16 mm	3/8 Circle Cutting	76 cm		FS 5001	FS 5002	
 16 mm	3/8 Circle round bodied	76 cm	FS 5081	FS 5082		
 16 mm	1/2 Circle Round Bodied	76 cm		FS 5049		
 20 mm	1/2 Circle Round Bodied	76 cm		FS 5086	FS 5087	
 25 mm	1/2 Circle Round Bodied	76 cm			FS 5070	
 25 mm	1/2 Circle Taper Cutting	76 cm				FS 5670
 26 mm	3/8 Circle Reverse Cutting	76 cm			FS 5028	FS 5029

NEEDLE LENGTH	NEEDLE DESCRIPTION	SUTURE LENGTH	SUTURE SIZE IN USP			
			2-0(3)	1-0(3.5)	1(4)	2(5)
 30 mm Heavy	1/2 Circle cutting	76 cm	FS 5072			
 30 mm	1/2 Circle cutting	76 cm		FS 5077		
 30 mm	1/2 Circle Round Bodied	76 cm	FS 5333	FS 5334		
		90 cm	FS 5331			
 36 mm	3/8 Circle cutting	76 cm	FS 5362	FS 5360	FS 5361	
 40 mm	1/2 Circle cutting	76 cm		FS 5064		
 40 mm	1/2 Circle reverse cutting	76 cm	FS 5067	FS 5068	FS 5069	
 40 mm	1/2 Circle Round Bodied	76 cm		FS 5075	FS 5063	FS 5073

NEEDLE LENGTH	NEEDLE DESCRIPTION	SUTURE LENGTH	SUTURE SIZE IN USP		
			2-0(3)	1-0(3.5)	1(4)
 45 mm	3/8 Circle Reverse Cutting	76 cm	FS 5036		
		90 cm		FS 5037	
 50 mm	1/2 Circle Reverse Cutting	76 cm	FS 5052		
 60 mm	3/8 Circle Cutting	76 cm			FS 5062

PRECUT SUTURES (NON NEEDED)

SUTURE LENGTH	PACKING	SUTURE SIZE IN USP			
		3-0(2)	2-0(3)	1-0(3.5)	1(4)
2 x 75 cm	Sterilised black braided silk in sterile foil over wrap. 2 x 75 cm box of 12's pack.	FS 212	FS 213	FS 214	FS 215

SILK REELS

SUTURE LENGTH 25 METERS	PACKING	SUTURE SIZE IN USP			
		3-0(2)	2-0(3)	1-0(3.5)	1(4)
6 Reels in a Box Non sterile	Surgical black braided 6 reels x 25 meters. (Non-sterile)	F 822	F 823	F 824	F 825